

Method and Nucleic Acids for the Detection of Microorganisms Relevant to Brewing

The invention relates to a method for the detection of microorganisms relevant to brewing, as well as to nucleic acids and combinations thereof which can be used in this method. The invention further relates to the use of the nucleic acids according to the invention or combinations thereof for the detection and/or for the identification and/or characterisation of different genera or species of microorganisms relevant to brewing.

Beer can be regarded as very stable microbiologically, and can only be spoilt by a relatively manageable number of bacteria. In order to discover contamination with these organisms as early as possible, an analytical system which allows rapid detection of the microorganisms in the matrix beer must be used, since countermeasures must be undertaken immediately.

The common feature of all microorganisms harmful to beer is the trace contamination of individual vessels (barrels, bottles) and their slow growth. In particular, microbiological culturing of the anaerobic microorganisms is very difficult. The beer-spoiling bacteria at present known are classed into the following genera: *Lactobacillus*, *Pediococcus*, *Pectinatus* and *Megasphaera*. Members of the *Selenomonas* and *Zymophilus* genera have not yet emerged as beer contaminants; however, contamination of beer and their subsequent growth in it cannot be ruled out.

The genus *Lactobacillus* describes Gram positive, non-sporulating, mostly immotile and chain-forming rods, which are long, thin and sometimes curved. Coccoid forms are also sometimes observed. Members of the genus *Lactobacillus* are microaerophilic, and some are anaerobic. They are cytochrome- and catalase-negative, their metabolism is fermentative and they require a complex nutrient medium. The molar G+C content of the DNA is between 32 and 53%.

As well as in beer, *Lactobacilli* are found in dairy and cereal products, in meat and fish products, in water, waste water, wine, fruit and fruit juices, acid-pickled vegetables, sauerkraut, silage and sourdough. Although they are a component of the normal oral, intestinal and vaginal flora of mammals, they are however seldom pathogenic (*Bergeys Manual of Syst. Microbiology*, 1984, p. 1209-1234). In beer, because of their metabolic

products, they lead to clouding and undesired flavour changes. Species relevant to beer spoilage are *Lactobacillus brevis*, *Lactobacillus lindneri*, *Lactobacillus casei*, *Lactobacillus paracasei*, *Lactobacillus coryniformis* and *Lactobacillus curvatus* (Back, Brauwelt, 1980, 120, p. 1562-1569).

The genus *Pediococcus* includes Gram positive, immotile and non-sporulating cocci. They form tetrads or occur as pairs. They are facultative anaerobes, and their oxygen sensitivity differs from species to species. *Pediococci* are cytochrome and catalase-negative and require a complex nutrient medium (*Bergeys Manual of Syst. Microbiology*, 1984, p. 1075-1079). They are used as starter cultures for the production of raw sausage products, they ferment various types of pickled vegetables and lead to the spoilage of foodstuffs (Firnhaber, Baumgart: *Mikrobiologische Untersuchung von Lebensmitteln*, 1993, p. 413-419, 115-117). The genus includes 8 species, and the species *Pediococcus damnosus* and *Pediococcus inopinatus* should be regarded as harmful to beer.

The genus *Pectinatus* includes the species *Pectinatus cerevisiiphilus*, *Pectinatus frisingiensis* and the strain *Pectinatus* sp. DSM 20764, not further taxonomically classified. All strains have been isolated from spoilt beer (Schleifer et al., *Int. J. of Syst. Bacteriology*, 1990, p. 19-27). These are slightly bent, non-sporulating rod-shaped bacteria. They have comb-like flagella, and are motile. They produce neither catalase nor cytochrome oxidase, and are obligate anaerobes. The molar G+C content is 38-41%. In the genus *Pectinatus*, and also in the genera *Megasphaera*, *Selenomonas* and *Zymophilus*, the cell wall is more similar to the Gram-positive bacteria than to the Gram-negative bacteria. Although the Gram staining is negative, they are taxonomically classified among the Gram-positive bacteria (Haikara, *The Prokaryotes*, 2nd Edition, Vol. II, 1991, p. 1993-2004).

The genus *Megasphaera* includes the species *Megasphaera elsdenii* and *Megasphaera cerevisiae*. Only *Megasphaera cerevisiae* is relevant to brewing, and is described as a Gram negative, strictly anaerobic, cytochrome- and catalase-negative, immotile and sometimes slightly stretched coccus, which occurs singly, in pairs or in short chains. The mean cell diameter is about 1.4 µm, and the molar G+C content 42.4-44.8%. Main metabolites are sulphur compounds, such as H₂S and volatile fatty acids. In beer,

contamination with *Megasphaera cerevisiae* leads to very marked changes in aroma and taste (Haikara, *The Prokaryotes*, 2nd Edition, Vol. II, 1991, p. 1993-2004).

Species of the genus *Selenomonas* are defined as obligate anaerobes, Gram negative, non-sporulating, slightly curved and motile rods. The molar G+C content is about 48-58% (Schleifer et al., *Int. J. of Syst. Bacteriology*, 1990, p. 19-27). *Selenomonads* are isolated from the stomach and intestinal tract and the dung of mammals. The genus includes 10 species (Hespell et al., *The Prokaryotes*, 2nd Edition, Vol. II, 1991, p. 2005-2013). Only *Selenomonas lacticifex* has been isolated from starter yeast, and is thus relevant to brewing. *Selenomonas lacticifex* has not yet emerged as a beer-spoiling bacterium; however, its growth in beer is possible, and hence it fulfils the definition of a beer-spoiling organism.

The species *Zymophilus paucivorans* and *raffinosisivorans* belong to the genus *Zymophilus* as Gram-negative, slightly bent, motile rods, which occur singly, in pairs or in short chains. The molar G+C content is about 38-41% . They are obligate anaerobes and have a fermentative metabolism. Both species are isolated from starter yeasts and brewery wastes; growth in beer has only been observed with *Zymophilus raffinosisivorans* (Schleifer et al., *Int. J. of Syst. Bacteriology*, 1990, p. 19-27).

On the basis of comparison of the 16S rRNA gene sequences, all the genera to be tested are classified among the Gram-positive bacteria with low G+C content. The genera *Pediococcus* and *Lactobacillus* are classified into the *Lactobacillaceae* family, and the genera *Pectinatus*, *Megasphaera*, *Selenomonas* and *Zymophilus* into the *Sporomusa* group. The *Sporomusa* group is also described as a group of the Gram-positive Eubacteriales with Gram-negative cell wall (Stackebrandt et al., *The Prokaryotes*, 2nd Edition, Vol. II, 1991, p. 25-26, 33).

A classical microbiological determination of the microorganisms described above can require up to 10 days. However, a markedly faster analysis is desirable, as otherwise unnecessary storage costs arise or the beer being tested has already been delivered. For these reasons, several rapid detection methods have already been developed. Thus, for example, organisms harmful to beer can be detected on the basis of their metabolic products (Haikara et al. *Microbiology*, 1995, 141, p. 1131-1137). Other indirect methods are turbidometry (Haikara et al., *ASBC*, 1990, p. 92-95) and measurement of

the ATP bioluminescence (Miller et al., *J. Inst. Brew.*, 1989, Vol. 95, p. 317-319). Detection by means of antibodies is also rapid and specific (Gares et al., *ASBC*, 1993, p. 158-163; Winnewisser et al., *Int. J. of Bacteriology*, 1995, 45, p. 403-405). With these methods, the disadvantage is that either non-specific parameters are tested or only one species or genus is detected in each case. Also, the equipment and staff cost is high. An overview of rapid methods for the detection of contaminants relevant to brewing is given by Dowhanick (*Cerevisia*, 1995, 20/4, p. 40-49).

The polymerase chain reaction (PCR; Mullis et al., see US 4,683,195, US 4,683,202 and US 4,965,188) is a rapid and effective method of specifically detecting organisms. A range of nucleic acids are known, through the use of which as primers and/or probes the specific detection of microorganisms relevant to brewing is possible. However, a disadvantage is that with the use of these nucleic acid molecules in an amplification or detection reaction, it is always only possible to detect a fraction of all possible microorganisms relevant to brewing. These PCR systems serve for the specific detection in each case only of individual species in an amplification reaction of the genera *Lactobacillus*, *Pediococcus*, *Pectinatus* and *Mega-sphaera* (Sakamoto US 5,869,642; Nietupski et al., US 5,705,339 and US 5,484,900; Tsuchia et al., JP 06141899A, JP 06113888A / *ASBC J.*, 1992, p. 64-67 / *ASBC J.*, 1993, p. 40-41; Yasui JP07289295A / *Can. J. Microbiol.*, 1997, 43, p. 157-163, Shimada et al., JP06090793; Alatossava et al. WO97/09448; Doyle et al., *J. of Ind. Microbiology*, 1995, 15, p. 67-70; DiMichele et al., *ASBC J.*, 1993, p. 63-66; Vogeser et al, *Brauwelt*, 1998, 24/25, p. 1060-1063). Further, the methods described for visualisation of the amplification products, such as, for example, agarose gel electrophoresis, present problems, as the carcinogenic and highly toxic ethidium bromide is used for staining the amplification products. These methods can only be automated with difficulty and the assessment of the agarose gels or the identification of the microorganisms on the basis of the length of the amplification products is sometimes not clear.

The problem to be solved by the present invention was, therefore, to provide a method and means which make possible a rapid test of beer and brewing raw materials for contamination with microorganisms, the test being required to detect the whole range of possible beer-contaminating microorganisms.

This problem is solved according to the invention by a process which comprises the following steps:

- (a) bringing the sample into contact with a combination of at least two first nucleic acid molecules (primers), which hybridise with a region of a microbial nucleic acid conserved in microorganisms relevant to brewing;
- (b) amplification of the microbial nucleic acid or a portion thereof to produce at least one amplification fragment;
- (c) bringing the amplification fragments obtained in step (b) into contact with at least one second nucleic acid molecule (probe), which specifically hybridises with at least one amplification fragment that comprises a sequence of the microbial nucleic acid specific for all microorganisms relevant to brewing or for one or several families, genera or species of microorganisms relevant to brewing; and
- (d) detection of at least one hybrid nucleic acid which consists of an amplification fragment and a second nucleic acid molecule introduced in step (c),

and by a nucleic acid molecule selected from:

- (i) a nucleic acid with a sequence according to SEQ ID NO 1-107 or a fragment thereof at least 10, preferably 15-30, nucleotides long;
- (ii) a nucleic acid which specifically hybridises with a nucleic acid according to (i);
- (iii) a nucleic acid which is at least 70%, preferably at least 90%, identical with a nucleic acid according to (i) or (ii), or
- (iv) a nucleic acid which is complementary to a nucleic acid according to (i) to (iii).

In the sequences according to SEQ ID NO 1-107, nucleotides are abbreviated as follows: G = guanosine, A = adenosine, T = thymidine, C = cytidine, U = uracil, i = inosine. In accordance with IUPAC, mixtures are abbreviated as follows: R = G or A, Y = C or T, K = G or T, W = A or T, S = C or G, M = A or C, B = C, G or T, D = A, G or T, H = A, C or T, V = A, C or G, and N = A, C, G or T.

For the determination of identity (in the sense of complete agreement, corresponding to 100% identity) with nucleic acid sequences according to (iii), partial sequences of a larger polynucleotide are considered. These partial sequences include 10 nucleotides and are identical when all 10 building blocks are identical in the two sequences compared. The nucleotides thymidine and uridine are to be regarded as identical. All possible fragments of a larger polynucleotide can be regarded as partial sequences.

Here 90% identity is present, when in the two sequences to be compared 9 out of 10 or 18 out of 20 nucleotides in one section are identical.

As an example, let us consider two polynucleotides which comprise 20 nucleotides and differ in the 5th element. In a sequence comparison, six 10-nucleotide ones are then found which are identical, and 5 which are not identical, as they differ in one element.

Otherwise, the identity can also be determined by degree, the unit being stated in percent. For determination of the degree of identity, partial sequences are also considered, which as a minimum include the length of the sequence actually used, e.g., as primer, or else 20 nucleotides.

As an example, polynucleotides A with a length of 100 nucleotides and B with a length of 200 nucleotides are compared. From polynucleotide B, a primer with a length of 14 nucleotides is derived. For the determination of the degree of identity, polynucleotide A is compared with the primer over its whole length. If the sequence of the primer occurs in polynucleotide A, but differs in one element, then there is a fragment with a degree of identity of $13/14 \rightarrow 92.3\%$.

In the second example, the whole of the aforesaid polynucleotides A and B are compared. In this case, all possible comparison windows of a length of 20 nucleotides are applied, and the degree of identity determined for them. Thus, if nucleotides 50-69 of polynucleotide A and B are identical with the exception of nucleotide No. 55, then for these fragments a degree of identity of $19/20 \rightarrow 95\%$ is found.

The method according to the invention can be carried out more rapidly than the previous microbiological detection methods, and makes it possible to detect several, preferably all, microorganisms relevant to brewing potentially present in a sample, such as, for example, even *Lactobacillus* species or members of the genera *Selenomonas* or

Zymophilus seldom arising as contaminants, for which hitherto no detection method existed. The detection is comprehensive and indicates all contamination risks in the brewery. By means of the method according to the invention, microorganisms relevant to brewing can be detected both in beer samples and also in raw material samples (barley malt, yeast, hops, water) or samples of intermediate products in beer production (e.g. mash, wort) even when the number of contaminating microorganisms is still low.

In this context, microorganisms relevant to brewing are understood primarily to mean bacteria and in particular the bacteria described above, *Lactobacillus brevis*, *Lactobacillus lindneri*, *Lactobacillus casei*, *Lactobacillus paracasei*, *Lactobacillus coryniformis*, *Lacto-bacillus curvatus*, *Pediococcus damnosus*, *Pediococcus inopinatus*, *Pectinatus cerevisii-philus*, *Pectinatus frisingiensis*, *Pectinatus* sp. DSM 20764, *Megasphaera cerevisiae*, *Selenomonas lacticifex*, *Zymophilus paucivorans* and *Zymophilus raffinivorans*, and also all microorganisms to be found in beer, which, while they do not belong to the aforesaid species, can nonetheless multiply in beer, for example, rare members of the *Lactobacillaceae* family, such as *Lactobacillus malefermentans*, *Lactobacillus buchneri*, *Lactobacillus parabuchneri*, *Lactobacillus sanfrancisco*, *Lactobacillus delbrueckii*, *Leuconostoc mesenteroides*, *Pediococcus pentosaceus* and *Lactococcus lactis*.

The microorganisms detectable by the method according to the invention are, thus, not limited to the microorganisms hitherto described as beer contaminants. Rather, the use of the nucleic acid molecules and the method according to the invention offers the possibility of recognising the presence of other microorganisms relevant to brewing, which have not previously been described as beer contaminants. A positive result at the level of higher taxonomic units (e.g. orders, families, genera) combined with a negative result at the level of the lower taxonomic units known to be relevant to brewing (e.g. species, subspecies, strains) indicates a contamination with such a non-typical microorganism relevant to brewing.

In a first step of the method according to the invention, the sample to be tested is brought into contact with a combination of at least two first nucleic acid molecules (primers). These nucleic acid molecules hybridise with a region of a microbial nucleic acid which is conserved in microorganisms relevant to brewing. The hybridisation takes place through pairing of the primer with regions of the microbial nucleic acid which have

an at least partly complementary base sequence. The term “conserved” characterises the evolutionary variability of nucleotide sequences for species of different taxonomic units. If corresponding sequence sections from at least two microorganisms relevant to brewing are compared, the sequence can be regarded as variable or as conserved. Comparison sequences which are at least 95% identical are described as conserved, and those which are less than 95% identical as variable. Thus, a region of a nucleic acid conserved in microorganisms relevant to brewing denotes a region which is at least 95% identical in all microorganisms relevant to brewing (as defined above).

In a preferred embodiment of the present invention, the conserved region occurs in a genome section which contains the bacterial 23S and 5S genes. This region includes the intergenic spacer between the genes for the 23S rRNA and the 5S rRNA and the bounding 23S and 5S rDNA genes, and includes both conserved sequence regions and also hypervariable (i.e., very organism-specific) sequence regions. Prokaryotic ribosomes as a rule contain three distinct nucleic acid components, which are generally known as 5S, 16S and 23S rRNA (ribosomal nucleic acid). The genetic information for these ribonucleic acids (rDNA) is typically arranged in the genome as a tandem. The typical organisation of such a unit is 16S-23S-5S, where the genes are connected to one another by short hypervariable intergenic regions, so-called spacers. The units are present several times in the genome, and the number of operons can vary from species to species. The high conservation of the DNA sequence in certain sections of the ribosomal DNA over the whole bacterial kingdom allows the design of non-specific oligonucleotides even without exact knowledge of the individual DNA sequences of the organisms to be investigated. The sequences according to SEQ ID NO 1-20 according to the invention (Table 1) are sequences of the 23S-5S intergenic spacer of microorganisms relevant to brewing, from which nucleic acid molecules for use in the method according to the invention can be derived.

The combination of at least two first nucleic acid molecules used in the first step of the method according to the invention is selected, such that they are usable as primers in an amplification reaction, i.e., one nucleic acid molecule hybridises onto a first conserved region of the first strand of the target DNA and the other nucleic acid onto a second conserved region of the DNA strand complementary to the first, wherein the desired target region of the DNA is included. Both nucleic acid molecules have a length of at

least 10 bp, preferably 15-30 bp. In a preferred embodiment of the invention, a combination of at least two nucleic acid molecules according to this invention is used. In a particularly preferred embodiment of the invention, a combination is used which includes at least one nucleic acid molecule with a sequence according to one of the SEQ ID NO 40 to 47 (Table 2) and at least one nucleic acid molecule with a sequence according to SEQ ID NO 48-54 or SEQ ID NO 55-59 or SEQ ID NO 60-72 (Table 2).

In a second step of the method according to the invention, the microbial nucleic acid or a portion thereof is amplified, whereby at least one amplification fragment is produced. Amplification is understood to mean the raising of the concentration of a nucleic acid or a portion thereof present in a reaction mixture. Processes used for the amplification of nucleic acids are for example the PCR (US 4,683,195, US 4,683,202 and US 4,965,188), the "self-sustained sequence replication" (EP 329,822), the "transcription-based amplification system" (EP 310,229) and the " β -RNA replicase system" (US 4,956,858). In a preferred embodiment of the present invention, the amplification comprises a polymerase chain reaction (PCR). In a further embodiment of the present invention, the amplification comprises a ligase-chain reaction or an isothermal nucleic acid amplification.

In a third step of the method according to the present invention, the amplification fragments obtained are brought into contact with at least one second nucleic acid molecule (probe). This nucleic acid molecule or these nucleic acid molecules hybridise specifically with at least one amplification fragment that comprises a sequence of the microbial nucleic acid which is specific for all microorganisms relevant to brewing or for one or several families, genera or species of microorganisms relevant to brewing, i.e., only occurs in members of these families or genera or in these species.

The double-strand formation of two identical or similar nucleotide fragments (DNA, RNA, PNA) is described as hybridisation. The term specific hybridisation is used when a stable hybrid nucleic acid between the oligonucleotide and the corresponding target DNA of the oligonucleotide exists, but not to other DNA than the target DNA. For the purposes of this invention, the feature "sequence which specifically hybridises with a sequence according to (i)" refers to a sequence, which under stringent conditions, hybridises with the sequence according to (i). For example, the hybridisations can be carried out at 50°C with a hybridisation solution consisting of 2.5 x SSC, 2 x Denhardt's solution, 10 mM

Tris, 1 mM EDTA pH 7.5. Suitable washing conditions are for example four times repeated 1-minute washings in 0.1 x SSC to 1.0 x SSC, 2 x Denhardt's, 10 mM Tris, 1 mM EDTA, pH 7.5 at 20-50°C.

In a preferred embodiment of the invention, one or several of the nucleic acid molecules according to the invention is used as a second nucleic acid molecule (probe).

Consensus probe is understood to mean a nucleic acid molecule which hybridises with highly conserved regions of a microbial nucleic acid and reacts with the amplification products of all microorganisms relevant to brewing. Nucleic acid molecules according to the invention which are usable as consensus probes have a sequence according to one of SEQ ID NO 40 to 72 (Table 2).

For the detection of a specific genus of microorganisms relevant to brewing, a nucleic acid molecule with a sequence according to one of SEQ ID NO 35 to 39 or SEQ ID NO 104 to 107 (Table 2) is preferably used. The genus specificity of a probe is defined as the ability of this probe to hybridise with the DNA of all isolates of as large as possible a group of members of the particular genus to be detected.

Species-specific nucleic acid probes are understood to mean nucleic acid molecules which hybridise with the DNA of all isolates of the particular species to be detected under the same stringency conditions. Species-specific nucleic acid molecules according to the invention with SEQ ID NO 21-22, SEQ ID NO 25-34, SEQ ID NO 73-78, SEQ ID NO 80-85 or SEQ ID NO 87-97 (Table 2) can be used.

The probes SEQ ID NO 23-24, SEQ ID NO 79, SEQ ID NO 86 and SEQ ID NO 98 to 103 are special cases. With the probes according to SEQ ID NO 23 and SEQ ID NO 79, strains of *Lactobacillus casei* and *Lactobacillus paracasei* ssp. *paracasei* can be detected. A probe according to SEQ ID NO 24 allows the detection of two subspecies of *Lactobacillus coryniformis* (*L. coryniformis* ssp. *coryniformis* and *L. coryniformis* ssp. *torquens*). With the probe SEQ ID NO 86, strains of the species *Pediococcus damnosus*, *Pediococcus inopinatus* and *Pediococcus parvulus* can be detected. With the use of these probes, other microorganisms relevant to brewing are not detected. Likewise, with the probes SEQ ID NO 98 to 103, all species of the *Lactobacillaceae* family relevant to brewing to be detected are detected, and other species and genera relevant to brewing are discriminated against.

In the last step of the method according to the invention, the detection of at least one hybrid nucleic acid which consists of an amplification fragment and a second nucleic acid molecule introduced in the preceding step takes place.

Preferably, first nucleic acid molecules (primers) and/or second nucleic acid molecules (probes) are at least 10 nucleotides, preferably 15-30 nucleotides long. In one embodiment of the present invention, the first and/or the second nucleic acid molecules are modified in that up to 20% of the nucleotides in 10 consecutive nucleotides, in particular 1 or 2 nucleotides of a block of 10 are replaced by nucleotides which do not occur naturally in bacteria.

The method according to the invention preferably includes the so-called consensus PCR. In this method, multiplication of the microbial nucleic acid or a portion thereof, and subsequent detection of these molecules by hybridisation with labelled specific probes take place. In the consensus PCR, nucleic acid molecules are used which make it possible to obtain an amplification product from several or, indeed, all of the relevant strains, subspecies, species or genera. The amplification does not lead to a differentiation of the microorganisms. The specificity of the detection is achieved through the subsequent hybridisation reaction with specific probes. In this way, microorganisms relevant to brewing can be simultaneously detected in a simple combination of amplification and detection reaction.

This kind of amplification and detection makes it possible to automate the detection reaction, so that a high sample throughput becomes possible. For example, a PCR-ELISA detection procedure can be used, in which the respective probes are bound in different wells of a microtitre plate, in which the hybridisation and the detection of the labelled amplification products then occurs. The detection can also be effected by the use of a microarray, on which several probes are immobilised, as a result of which the detection reaction can be carried out quickly and at no great cost.

In a preferred embodiment of the invention, the second nucleic acid molecule (probe) is modified or labelled in such a way that it can produce a detectable signal. The modification or labelling is selected from (i) radioactive groups, (ii) coloured groups, (iii) fluorescent groups, (iv) groups for immobilisation on a solid phase and (v) groups

which permit an indirect or direct reaction, especially with the aid of antibodies, antigens, enzymes and/or substrates with affinity to enzymes or enzyme complexes.

For the purposes of this invention, labelling indicates directly or indirectly detectable groups or groups for immobilisation on a solid phase, which are attached to the nucleic acid molecule. Directly detectable are metal atoms, radioactive, coloured or fluorescent groups. Indirectly detectable are immunologically or enzymatically detectable groups, for example, antigens and antibodies, haptens or enzymes or enzymatically active parts of enzymes. These indirect groups are detected in subsequent reactions. Preferred are haptens which are coupled to an oligonucleotide and which are detected in a subsequent antibody reaction.

The nucleic acid molecules according to the invention can be used for the detection and/or for the identification and/or characterisation of bacteria relevant to brewing. The primers and/or probes described herein can also be used in the detection of the described microorganisms in drinks other than beer, in other samples from the brewing sector, such as for example in raw materials, starter yeast, environmental samples, in other foodstuff samples or in clinical samples, etc.

Examples:

Example 1: Determination of the DNA target sequence of the bacteria harmful to beer and closely related species

By sequence comparison of known 23S rDNA and 5S rDNA sequences (GenBank Sequence Database of the National Center of Biotechnology Information: NCBI), conserved gene regions were identified, which serve as hybridisation sites for the primers used for the sequencing. From pure cultures of the bacteria listed in Table 1, genomic DNA was isolated by known standard methods. With primers which hybridise in highly conserved regions, amplification products of all bacteria to be detected were obtained in a PCR. The following primers were used for the amplification and the subsequent sequencing:

Primer 1 = SEQ ID NO 47:

5'-AAG TGC TGA AAG CAT CTA AG-3'

Primer 2 = SEQ ID NO 55:

5'-GGC RRY GTC TAY TYT CSC-3'

Composition of the PCR:

Genomic DNA (10 – 100 ng)	1.00 μ l	
H ₂ O	16.85 μ l	
Buffer (10 x)	2.50 μ l	1 x
dNTP (10 mM)	0.50 μ l	200 μ M
Primer 1 = Seq ID NO 48 (5 μ M)	1.50 μ l	0.30 μ M
Primer 2 = Seq ID NO 49 (5 μ M)	1.50 μ l	0.30 μ M
MgCl ₂ (50 mM)	1.00 μ l	2.00 mM
Taq-polymerase (5 U/ μ l)	0.15 μ l	0.03 U/ μ l
Σ	25.00 μ l	

Temperature profile:

5 mins	95°C	
30 secs	95°C	
30 secs	50°C	x 38
30 secs	72°C	
5 mins	72°C	

These amplification products were purified via an agarose gel and by a subsequent treatment with the QIAquick PCR Gel Extraction Kit (Quiagen Co.) and sequenced in the Long Read Sequencer Model 4000L (LI-COR Co.) with the aforesaid primers, which are provided with an IRD-800 label. The resulting sequences of the 23S/5S rDNA spacer regions of the bacteria relevant to brewing and the phylogenetically closely related species were compared with one another and sequence regions identified which:

- 1.) are to be found in all species of the particular genus to be detected and at the same time differ from those of other genera or species,
- 2.) are only to be found in the particular species to be detected, but differ from other bacteria to be detected and not to be detected.

In the sequence regions described under 1.), hybridisation sites of genus-specific oligonucleotides were defined, and in the sequence regions described under 2.), the binding sites of species-specific oligonucleotides were defined.

Example 2: Detection of Bacteria Harmful to Beer by the Polymerase Chain Reaction

I. Amplification

Genomic DNA was isolated from pure cultures of the bacteria listed in Table 1 by known standard methods. Decimal dilutions from 1 fg/μl to 1 pg/μl of these preparations were then used in a PCR with the following composition:

Primer 3 = SEQ ID NO 46:

5'-AAG GGC CAT CRC TCA ACG G -3'

Primer 4 = SEQ ID NO 48:

5'-TGT GTT CGT CAT GGG AAC AGG TG -3'

Genomic DNA	1.00 μl	4.00 μl	
H ₂ O	16.60 μl	66.40 μl	
Buffer (10 x)	2.50 μl	10.00 μl	1 x
dNTP (10 mM)	0.50 μl	2.00 μl	0.20 mM
Primer 3 = Seq ID NO 21 (5 μM)	1.50 μl	6.00 μl	0.30 mM
Primer 4 = Seq ID NO 22 (5 μM)	1.50 μl	6.00 μl	0.30 mM
digoxigenin labelled			
DMSO (100%)	0.25 μl	1.00 μl	1.00 %
MgCl ₂ (50 mM)	1.00 μl	4.00 μl	2.00 mM
Taq-polymerase (5 U/μl)	0.15 μl	0.60 μl	0.03 U/μl
Σ	25.00 μl	100.00 μl	

The PCR was performed under the following conditions in the Mastercycler® (Eppendorf Co.) according to the following temperature profile:

5 mins	95°C	x 38
30 secs	95°C	
45 secs	55°C	
90 secs	72°C	
5 mins	72°C	

Primer 3 (SEQ ID NO 46) was determined by sequence comparison of known 23S rDNA sequences (GenBank Sequence Database of NCBI). It hybridises onto highly conserved sequence sections in the 23S rDNA gene region. The binding site lies outside the region sequenced with the primers SEQ ID NO 48 and 49.

Primer 4 (SEQ ID NO 48) was determined on the basis of our own sequence data. The hybridisation site of primer 2 lies adjacent to the intergenic 23S/5S spacer in the 5S rDNA region.

II. Detection by PCR-ELISA

The detection is effected by PCR-ELISA. For this, per probe used, 5 µl of amplification product are treated with 5 µl of denaturation buffer (125 mM NaOH, 20 mM EDTA, pH 14) and incubated for 15 mins at room temperature. Each time, 2 pmoles of the particular biotinylated probe are pipetted into 100 µl of hybridisation buffer (2.5 x SSC, 2 x Denhardt's solution, 10 mM Tris, 1 mM EDTA, pH 7.5) and transferred to the wells of a microtitre plate coated with streptavidin and preincubated at the hybridisation temperature of 50°C. After the denaturation, the denaturation mixture is pipetted into the hybridisation mixture. Next the mixture is incubated for 30 minutes at hybridisation temperature. If the hybridisation is complete, the hybridisation mixture is removed and the plate washed 4x with 200 µl of wash buffer 1 (WB1: 0.1 x SSC, 2 x Denhardt's, 10 mM Tris, 1 mM EDTA, pH 7.6) for 1 min. each time at hybridisation temperature. Next, 100 µl of a solution of a horseradish peroxidase conjugated anti-digoxigenin antibody diluted according to the manufacturer's instructions is added (Boehringer Mannheim). The conjugate is diluted in wash buffer 2 (WB2: 100 mM Tris, 150 mM NaCl, 0.05% Tween 20, 0.5% blocking reagent, 100 µg/ml herring sperm, pH 7.6). Next, the antibody incubation is performed at 37°C for 30 mins. After this, the plate is washed four times with 200 µl of WB2 (at room temperature). After the washing, 100 µl of POD

substrate (Boehringer Mannheim) are added and the mixture incubated for 20 mins at RT. Next the colour reaction is stopped with 100 µl of 0.5M H₂SO₄ and estimated at 450 nm.

III. Assessment

According to the detection protocol described above, the detection was performed for all bacteria and bacteria groups investigated, using the corresponding genus- and species-specific probes. Genus-specific probes are SEQ ID NO 35 for *Pediococcus*, SEQ ID NO 36 for *Pectinatus*, SEQ ID NO 37 for *Megasphaera*, SEQ ID NO 38 for *Selenomonas* and SEQ ID NO 39 for *Zymophilus*. Species-specific probes are SEQ ID NO 21 for *Lactobacillus brevis*, SEQ ID NO 22 for *Lactobacillus lindneri*, SEQ ID NO 23 for *Lactobacillus casei* + *paracasei*, SEQ ID NO 24 for *Lactobacillus coryniformis*, SEQ ID NO 25 for *Lactobacillus curvatus*, SEQ ID NO 26 for *Pediococcus damnosus*, SEQ ID NO 27 for *Pediococcus inopinatus*, SEQ ID NO 28 for *Pectinatus cervisiophilus*, SEQ ID NO 29 for *Pectinatus frisingiensis*, SEQ ID NO 30 for *Pectinatus* sp. DSM20764, SEQ ID NO 31 for *Megasphaera cerevisiae*, SEQ ID NO 32 for *Selenomonas lacticifex*, SEQ ID NO 33 for *Zymophilus paucivorans* and SEQ ID NO 34 for *Zymophilus raffinovorans*.

As controls, the consensus probes SEQ ID NO 40 and 41 were used, which hybridise with the amplification products of all the species to be detected. Further possible binding sites for consensus probes are SEQ ID NO 42-45. The probes of SEQ ID NO 40 to 45 were determined by sequence comparison of known 23S rDNA and 5S rDNA sequences (GenBank Sequence Database, NCBI).

If the extinction measured for a 1 fg quantity of genomic DNA used in the PCR was greater than 1, the result was assessed as positive. The results of the PCR-ELISA are presented in Table 3.

Table 1

SEQ ID NO	Source			Description	Sequence
	Genus	Species	Strain		
1	Lactobacillus	brevis	DSM 20054	23S-spacer-5S	5'-TATATGGAAG TAAGACCCCT GAGAGATGAT CAGGTAGATA GGCTGGAAGT AGCAGGCGCG TGAGGGGTGG AGCGGACCAG TACTAATCGG TCGAGGACTT AACCAAGTCA ACAACGTAGT TGTTTCGAGA ATAATTGAAT AATATCTAGT TTTGAGGGAA GAAGTCTCT TATAGTGTGG TGGCGATAGC CTGAAGGATA CACCTGTTCC CATGCCGAAC ACAGAAGTTA AGCTTCAGCA CGCCGATAGT AGTTGGGGGA TCGCCCC-3'
2	Lactobacillus	lindneri	DSM 20690	23S-spacer-5S	5'-CCATTCTTAT ATGGAAGTAA GACTCCTGAA AGATGATCAG GTCGATAGGT TAGAAGTGGA AGCATAGTGA TATGTGAAGC GGACTAATAC TAATCAGTCG AGGACTTAAC CAAGGAAGAC ACAGGTTAA ATCAAAGTTG AACAGAGAAG ATATTACTA GTTTTGAGAG AACGAAGTTC GCTCAGGCTT ATGAAAAATA AGCATAGTGT GTTGGCGATA GCCTGAAGGA TACACCTGTT CCATGCCCCG ACACAGAAGT TAAGCTTCAG CACGCCAAA GTAGTTGGGG GATCGCCCCC TCGCAGGATA GGACGATGGT CATAGC-3'
3	Lactobacillus	casei	DSM 20011	23S-spacer-5S	5'-CCATTCTTAT ATGGAAGTAA GACCCCTGAG AGATGATCAG GTAGATAGGC TGGAAGTGGA AGTGCAGCGA TGCAATGGAGC GGACCAGTAC TAATCGGTCTG AGGACTTAAC CAAGTAGAGC GTGAGCAGGA GCGCTTAGAA ACCGAGGCAT AAGCGGCGCT GAGTTCGTTG GCCGGTTTT GGCCAATGGA TTCAGGTTTC TTATGTGGAG GTTCTGCGA CTGCGAACGC GTTTCGATGA AATACACTGG TTCCCGACAA CACAAAAACA ACAATGATAG CCAGTTTGA GAGCGCAAGG TTCTCATAAG TGTGTTGGCG ATAGCAAGAA GGATACACCT GTTCCCATGC C-3'
4	Lactobacillus	paracasei ssp. paracasei	DSM 20008	23S-spacer-5S operon 1	5'-CCATTCTTAT ATGGAAGTAA GACCCCTGAG AGATGATCAG GTAGATAGGC TGGAAGTGGA AGTGCAGCGA TGCAATGGAGC GGACCAGTAC TAATCGGTCTG AGGACTTAAC CAAGTAAGAG TGTGAGCAGG AGCGGTTAGA AACCGGAGCA TAAGCGGGCC TGAGCGTGAT GGCCGGGCTT TGGCCATTGC GGTACGGTTC CTTATGTGCA GGTTCCTGCG ACTGCGAACA CGTTTCGATG ACAAGTACGT TAAGTTCAAG GCAGCAATTA AACATGATA GCTAGTTTTC AGAGCGCAAA GTTCTCATAA GTGTGTTGGC GATAGCAAGA AGGATACACC TGTTCCCATG CCGAACACAG AAGTTAAGCT TCTTCACGCC GAGAGTAGTT GGTGGGAAAC TGCCTGCGAG GATA-3'

SEQ ID NO	Source		Description	Sequence	
	Lactobacillus	paracasei ssp. paracasei		DSM 20008	23S-spacer-5S operon 2
5	Lactobacillus	paracasei ssp. paracasei	DSM 20008	23S-spacer-5S operon 2	5'- CCATTCCTAT ATGGAAGTAA GACCCCTGAG AGATGATCAG GTAGATAGGC TGGAAAGTGA AGTGAAGCGA TGCAATGGAG GGACCAGTAC TAATCGGTGC AGGACTTAAC CAAGTAAGCG TGCAAGCAGG AGCAGGTTTC TGCAGACTGC AACACATTC GATGACAACT ACCTTAAGTT CAAGGCAGCA ATTAACAGAT CATAGCCAGT TTTGAGAGCG CAAAGTTCTC ATAAGTGTGG TGGCGATAGC AAGAAGGATA CACCTGTTC CATGCCGAAC ACAGAAAGTTA AGCTTCTTCA CGCCGAGAGT AGTTGGTGGG AAACGTCCCTG CGAGGATA-3'
6	Lactobacillus	coryniformis ssp. coryniformis	DSM 20001	23S-spacer-5S	5'- CTCGAGTTGA GATTCCCAT TCCTTTATGG AAGTAAGACC CCTGAGAGAT GATCAGGTAG ATAGTTTGA AGTGGACGTG CCGTGAGGCA TGGAGCGGAC CAATACATA CCGTCGAGGA CTTAACCAAG TAGCATGTAC GTAGTGTAG TTTAAGGGCA AAGAAATGAA TATCCAGTTT TGAGAGCGCA ACGTTCTCAG AAAGTGGTGT GGTGGCGATA GCAAGAAGGA TACACCTGTT CCCATGTCCA ACACAGAAGT TAAGCTTCTT AGCGCCGAGA GTAGTTGGGG GAGCACCCCC TGGGAGGATA GGACGAT-3'
7	Lactobacillus	coryniformis ssp. torquens	DSM 20004	23S-spacer-5S	5'- CTCGAGATGA GATTCCCAT TCCTTTATGG AAGTAAGACC CCTGAGAGAT GATCAGGTAG ATAGTTTGA AGTGGACGTG CCGTGAGGCA TGGAGCGGAC CAATACATA CCGTCGAGGA CTTAACCAAG TAGCATGTAC GTAGTGTAG TTTAAGGGCA AAGAAATGAA TATCCAGTTT TGAGAGCGCA ACGTTCTCAG AAAGTGGTGT GGTGGCGATA GCAAGAAGGA TACACCTGTT CCCATGTCCA ACACAGAAGT TAAGCTTCTT AGCGCCGAGA GTAGTTGGGG GAGCACCCCC TGGGAGGATA GGACGAT-3'
8	Lactobacillus	curvatus	DSM 20019	23S-spacer-5S	5'- ACGCCTCGAG ATGAGATTTC CCATTCCTTT ATGGAAGTAA GACCCCTGAA AGATGATCAG GTAGATAGGC TAGGAGTGA ACTACAGCGA TGATGAGGC GGACTAGTAC TAATCGGTGC AGGACTTAAC CAAAGGTGCA ATGTTAGGCT TTTGAAATGA AATATTACTT ATTATGCAGT TTTGAGAGCA GAAGTCTTT CTCAGTGGC AAGCACAAA TAGTGTGGTG GCCATAGCAA CCAAGATACA CCTGTTCCCA TGTCGAACAC AGAAGTTAAG CTTCTTAGCG CCGATAGTAG TTGGTGGGAA ACTACCTGCG AGGATAGGAC GATGAT-3'
9	Pediococcus	damnosus	DSM 20331	23S-spacer-5S	5'- GATGAGATTT CCCATTCAT TTAATGGAAGT AAGACCCCTG AGAGATGATC AGGTAGATAG GTTGGGAGTG GAAGTGTAGT GATACATGGA CGGACCAAT ACTAATCGGT CGAGGACTTA ACCACAAAAGT GGTGTTCTCA AGAAGAGAT TCGATATTAT TTAGTTTGA GAGAATAAAT TTTCTTTCACA CGAGCCGCGT AAGTGGATCG GAGAAGTGTG GTGACGATAT TGAGAAGGAT ACACCTGTTT CCATGTCGAA CACAGAAGTT AAGCTTCTTA ACGCCGAGAG TAGTTGGGGG ATCGTCCCT GCGAGGATAG GACGATGCTC AATAG-3'

SEQ ID NO	Source		Description	Sequence	
	Pediococcus	inopinatus		DSM 20285	5'- 23S-spacer-5S
10				5'- 23S-spacer-5S	AGATGAGATT TCCCATTTCCA TTTATGGAAG TAAGACCCCT GAGAGATGAT CAGGTAGATA GGTGGGAGT GGAAGTGTAG TGATACATGG AGCGGACCAA TACTAATCGG TCGAGGACTT AACACACAAAG TGGTGTCTC AAAGAGAAGA TTTTCGATATT ATTTAGTTTT GAGAGAATAA ATTTCTTTCA CACGAGCCGC GGAAGTGGAT CGGAGAAGTG TGGTGACGAT ACTGAGAAGG ATACACCTGT TCCCATGTGC AACACAGAAG TTAAGCTTCT TAACGCCGAG AGTAGTTGGG GGATCGCTCC CTGCGAGGAT AGGACG-3'
11	Pectinatus	cerevisiiphilus	DSM 20467	5'- 23S-spacer-5S	AAGTGCTGAA AGCATCTAAG CGTGAACCT GCCTTAAGAT GAGGTTTCCC AGAGCCGTAA GGTGTGAAG GCACCTTGAA TAAGACGAGG TAGATAGGCC GGGAGTAGAA GTACAGTAAT GTACGAAGCG GACTGTACT AATAAGCCGA GAGCTTAAC TAAATCATC GAAAAAATG TTTGGTCTGA GATTTCTTCT GTGAAGTTTT GAGTGTCAA GACACTCTGG TTGAAGGGA GGGAACTGA GAGCGTAAAA CTGCGGACTT TGGCTCAAA AGTTAAAGCA TCTGGTGAGG ATACTGGAT GGATCCACCT GTTCCCATTC CGAACACAGT AGTTAAGCAT CCACAGGCTG AAGGTACTTG GGGGGCGACC CCTTGGGAAA ATAGGACACT GCC-3'
12	pectinatus	frisingensis	DSM 6306	5'- 23S-spacer-5S	AAGTGCTGAA AGCATCTAAG CGTGAACCA GCTTTAAGAT GAGGTTTCCC AGAACGCAAG TTTGGAAGC ACCTTGAAG AGACGAGGTA GATAGCCCGG GAGTGAAGT ATGGTGACAT ATGAAGCCGA CTGGTACTAA TAAGCCGAGA GCTTAAC TTG ATTTTCATCAA AAAAGAGAAA TGTTTGGTCA GAGATTTTCT TCTGTGAAGT TTTGAGTGTG CAAGAACACT CGAGGTATA TAGGTAAGG AAAAGCAGCA GATAAGTTTC TATATACCCG TATATACCCG CTGAGGTGCT GAGGCACTGA AGGCCAGAAC ATCTGGTGCC GATACCTGGA TGGATCCACC TGTTCCTCAT CCGAACACAG TAGTTAAGCA TCCACAGGCC GAAGGTACTT GGGGGGCAGC CCCCTGCGAA AATAGGACAC CGCC-3'

SEQ ID NO	Source		Description	Sequence
13	Pectinatus	sp.	DSM 20764 23S-spacer-5S operon 1	50 AAGTGCTGAA AGCATCTAAG CGTGAAACCT GCCTTAAGAT GAGGTTTCCC 100 AGAGCCGTAA GGCTTGGAG GCACCTTGAA GATGACGAGG TAGATAGGCC 150 GGGAGTAGAA GTATGGTGAC ATACGAAGCG GACTGGTACT AATAAGCCGA 200 GAGCTTAACT TAAATTCATC TATAAATGTT TGGTCCCTGAT TTCTTCTGTG 250 AAGTTTGGAG TGTGCAAGAT CACTCATGAA AGTATATAGG TAAAGGGAAA 300 GCAGCAGATT AGTTCCCTGGT TTACTTTTATA TATGAGCACT AAGGTGCAGA 350 AAAGAACGTT TGAGGAAACG CGCGCTTCGT AAATCACTT TGCCTGCTGA 400 TTATCTCAAT GCTAAAGCAT TAAGATAATT TTAGAGGAAA CGCGCTTCA 450 CTAGCGTTCA CTCTCGGTAC TTTATTTCTA AGTGCTGAAG CACTAAGAAG 500 GGCAAGGAAA CGCGTCGTTT GCGATGCTCA CTTTGGGTAC TTCTATCTTA 550 GACTGCTAAA GCAGTAAGAT CTGAAGCATC TGGTGGCGAT ACCTGGATGG 600 ATCCACCTGT TCCCATTCG AACACAGTAG TTAAGCATCC ACAGGCCGAA GGTACTTGGG GGCAGCCCC CTGCGAGAGT AGGACATCGC C-3'
14	Pectinatus	sp.	DSM 20764 23S-spacer-5S operon 2	50 AATGCTGAA AGCATCTAAG CGTGAAACCT GCCTTAAGAT GAGGTTTCCC 100 AGAGCCGTAA GGCTTGGAG GCACCTTGAA GATGACGAGG TAGATAGGCC 150 GGGAGTAGAA GTATGGTGAC ATACGAAGCG GACTGGTACT AATAAGCCGA 200 GAGCTTAACT TAAATTCATC TATAAATGTT TGGTCCCTGAT TTCTTCTGTG 250 AAGTTTGGAG TGTGCAAGAT CACTCATGAA AGTATATAGG TAAAGGGAAA 300 GCAGATTAGT TCCTTGGTTA CTTTATATAT GAGCACTAAG GTGCAGAAAA 350 GAACGCTAA GGAACCGCGG CGTTCGTAGG CTCACCTCTG GTACTTCTATC 400 TCTAGACTGC TAAAGCAGTA AGATCTGAAG CATCTGGTGG CGATACCTGG 450 ATGGATCCAC CTGTTCCCAT TCCGAACACA GTAGTTAAGC ATCCACAGGC CGAAGGTACT TGGGGGGCAG CCCCCTGCGA AAGTAGGACA CCGCC-3'
15	Megasphaera	cerevisiae	DSM 20462 23S-spacer-5S operon 1	50 GCATCTAAGC GTGAAACCAG CCTAGAGATG AGGTTTCTCA TTACGAAAGT 100 AAGTAAGTTC CCATGAAGAC GACATGGTAG ATAGGCCGGG AGTGGACGTA 150 CAGTAATGTA TGGAGCGGAC CGGTACTAAT AGACCGAGGA CTTGACTTAA 200 GCAGGGAACC CATTTTAAAG AAGCGAAGCG ACGCATAAAA TGGAGTGAAT 250 CGCTTATACC GAATCGCAGA TTTCGGTAAAG CAGCGAGAGAA TACCAATGCA 300 GCGGCAACAC CAGTTAGCAT AAATAAGCGG GATTCGGAGT GGTGAGGGA 350 GTTTCGTAGC AGCGTAGGCT AACCAACCA CCGCTTTCGA AGAAGGCGAA 400 TGGTTTGAAA AAGAGTACAT GCGAAGAAAC GACGAAAGAC TCACAACCAA 450 AACATACAAA CTAAGTAGAT GACATTAGAG TCACACCGAT TGTTAAGATC 500 CGAAATACTT TTCGATGTAG TTGTCAGGAT ACGAATCCTG AAACGAATTC AGTGGTGATG GCTGCAGGGA TCCACCTGTT CCCATACCGA ACACAG-3'

SEQ ID NO	Source		Description	Sequence
16	Megasphaera	cerevisiae	DSM 20462 23S-spacer-5S operon 2	5'-GCATCTAACC GTGAACCACG CCTAGAGATG AGGTTTCTCA TTACGAAAGT AAGTAAGGTC CCATGAAGAC GACATGGTAG ATAGGCCGGG AGTGACGTA CAGTAATGTA TGGAGCGGAC CGGTACTAAT AGACCGAGGA CTTGACTTAA GCAAAGAAGC AATAGAAAAG ACCATGTAGA TGGTGTAAAG GTTAGACGGG TAGTTAAGGT CCGAAATACT TTTCGATGTA GTTGTCAGGA TACGAATCCT GAAACGAATT CAGTGGTGAT GGCTGCAGGG ACCACTGTT CCCATACCGA ACACAG-3'
17	Selenomonas	lactificifex	DSM 20757 23S-spacer-5S operon 1	5'-AAGTGCTGAA AGCATCTAGG CGTGAAGCCT GTCCCGAGAT GAAGTATCTC ATGGAGTAAT CCAGTAAGAT TCCTTGAAGA AGACAAGGTA GATAGTTGG GAGTGAAGC ATCGTAAGGT GTTCAGCGGA CCAACTACTAA TAAATCGAGG GCTTAACCTT ACAGACCTGT CCAAGAGCG AGCGGATCG GGTAAACAGG CGTATGCGAA AACATCCAA GAATCGATC CAAAGGGCG AGATGATTGG CAGATGTTGA CCGCTAATAA TCTAGAATGT TTCGATACAA TTTTCTTCT GTATAGTTT GAGTGGACAT CGTTCAATCA ATAATATCCA GTACGATAG CTGAGTGGTA CCACCTGTTT CCATACCGAA CACAGTAGT AAGCACTCAT ACGCCGAAAG TACTTGTCTG GAAACGGGCT GCGAGAATAG GACGTCGCC -3'
18	Selenomonas	lactificifex	DSM 20757 23S-spacer-5S operon 2	5'-AAGTGCTGAA AGCATCTAAG CGTGAAGCCT GTCCCGAGAT GAAGTATCTC ATGGAGTAAT CCAGTAAGAT TCCTTGAAGA AGACAAGGTA GATAGTTGG GAGTGAAGC ATCGTAAGGT GTTCAGCGGA CCAACTACTAA TAAATCGAGG GCTTATCTTA ATAATCTAGA ATGTTTCGAT ACAATTTTC TTCTGTATAG TTTTGAGTGG ACATGGTTCA TTCAATAATA TCCAGTGACG ATAGCTGAT GGTACCACCT GTTCCCATAC CGAACACAGT AGTTAAGCAC TCATACGCCG AAAGTACTTG TCTGAAACG GGCTGCGAAA ATAGGACGCC GCC-3'
19	Zymophilus	raffinivorans	DSM 20765 23S-spacer-5S	5'-AAGTGCTGAA AGCATCTAAG CGTGAAGCCT GCCTTAAGAT GAGGTTTCTC ACAGAGCAAT CTGGTAAGAC CCCTTGAAGA AGACAAGGTA GATAGTTCGG GAGTGAAGC GCAGTAATGT GTGCAGCGGA CCGATACTAA TAGGTCGAGG GCTTGACTTA AGCCAGAAC GAAAACATAA ATGCGAACAT TTCTTCTTC TGATAGTTT TGAGAGAAC AACTCTTAAG ATGGAGTAG CTGAGGCGAA AGCGGAAGGC AGCGATATCT AAAAAAGAA TATCTGGTAG TGATAGCCAA GTGGACCCAC CTGTTCCCAT ACCGAACACA GTAGTTAAGC ACTGAACGT CGAAAGTACT TGGGTGAAA CGCCCTGCGA AAATAGGACA CCGCC-3'

SEQ ID NO	Source		Description	Sequence	
	Zymophilus	paucivorans			
20			23S-spacer-5S	5'- AAGTGCTGAA AGCATCTAAG CGTGAAACCA GCCTTAAGAT GAGTTTCTC ACAGAGCAAT CTGGTAAGAC CCCTTGAAGA AGACAAGGTA GATAGGTCGG GAGTGAAGC GCAGTAATGT GTGTAGCCGA CCGATACTAA TAGGTCGAGG GCTTGACTTA AAGCCAGAAC GAATTTCTAA ATGCGAACAT TTCTTTCTTC TGTATAGTTT TGAGAGAACA GACTCTTAAG ATGAGCAGTC TGAGGCGAAA GCTAAAGGCA GCGATATCTA AAAAAAGAA TATCTGGTAG TGATAGCCAA GTGGACCCAC CTGTTCCCAT ACCGAACACA GTAGTTAAGC ACTTGAACGT CGAAAGTACT TGGGTGGAAA CGCCCTGGGA AAATAGGACA CCGCC-3'	50 100 150 200 250 300 350

Table 2

SEQ ID NO	Description		Sequence	
21	Lactobacillus brevis	Specific probe	5'-	CCAAGTCAACAACGTAGTTGT
22	Lactobacillus lindneri	Specific probe	5'-	GACACAGGGTTAAATCAAAGTTG
23	Lactobacillus casei Lactobacillus paracasei ssp. paracasei	Specific probe	5'-	AGGTTTCTGCGACTGCGAAC
24	Lactobacillus coryniformis ssp. coryniformis Lactobacillus coryniformis ssp. torquens	Specific probe	5'-	ATGTACGTAGTGTTAGTTAAAGGCG
25	Lactobacillus curvatus	Specific probe	5'-	CTTCTCAGTGCACAAGCACA
26	Pediococcus damnosus	Specific probe	5'-	GTGTTCTCAAGAGAAGGATTCTG
27	Pediococcus inopinatus	Specific probe	5'-	GTTCTCAAAGAGAAGATTTCGATATTA
28	Pectinatus cerevisiiphilus	Specific probe	5'-	TGAGAGCGTAAAACTGCGACTT
29	Pectinatus frisingensis	Specific probe	5'-	CAGATAAGTTTCTCTGGTTACTG
30	Pectinatus sp. DSM 20764	Specific probe	5'-	CACTAAGGTGCAGAAAAAGAACGT
31	Megasphaera cerevisiae	Specific probe	5'-	CTTTTCGATGTAGTTGTCAGGATACG
32	Selenomonas lacticifex	Specific probe	5'-	GTTTCATTCAATAATATCCAGTGACG
33	Zymophilus raffinosisivorans	Specific probe	5'-	AACTCTTAAGATGGAGYAGTCTG
34	Zymophilus paucivorans	Specific probe	5'-	ACTCTTAAGATGAGCAGTCTGA
35	Pediococcus	genus-specific probe	5'-	AGTSTAGTGATACATGGAGCG
36	Pectinatus	genus-specific probe	5'-	GTGAAGTTTGTAGTTGCAAGA
37	Megasphaera	genus-specific probe	5'-	GACCGAGGACTTGACTTAAGCA
38	Selenomonas	genus-specific probe	5'-	TCCAGTGAACGATAGCTGAGT
39	Zymophilus	genus-specific probe	5'-	AAGAATATCTGGTAGTGATAGCCAA

Table 2 (Cont.)

SEQ ID NO	Description		Sequence	
40	consensus sequence	5'-	GTCGTGAGACAGTTCGGTC	-3'
41	consensus sequence	5'-	CYTAGTACGAGAGCCGRR	-3'
42	consensus sequence	5'-	GCTACCTGGGATAACAGGC	-3'
43	consensus sequence	5'-	ATCGACGGGAGGTTTSSCAC	-3'
44	consensus sequence	5'-	CACCTCGATGTCGGCTCRTC	-3'
45	consensus sequence	5'-	CCAAGGTTGGGCTGTTTC	-3'
46	consensus sequence	5'-	AAGGCCATCRCTCAACGG	-3'
47	consensus sequence	5'-	AAGTGTGAAGCATCTAAG	-3'
48	consensus sequence	5'-	TGTGTTGCGiATGGGAACAGGTG	-3'
49	consensus sequence	5'-	TGTGTTGCGAATGGGAACAGGTG	-3'
50	consensus sequence	5'-	TGTGTTGCGAAATGGGAACAGGTG	-3'
51	consensus sequence	5'-	TGTGTTGCGATATGGGAACAGGTG	-3'
52	consensus sequence	5'-	TGTGTTGCGCATGGGAACAGGTG	-3'
53	consensus sequence	5'-	TGTGTTGACATGGGAACAGGTG	-3'
54	consensus sequence	5'-	GGCRRYGTCTCTAATYTCSC	-3'
55	consensus sequence	5'-	GGCAGTGTCTACTTTTCCC	-3'
56	consensus sequence	5'-	GGCAGCGTCTACTTTTCGC	-3'
57	consensus sequence	5'-	GGCAGTGTCTACTTTTCGC	-3'
58	consensus sequence	5'-	GGCAGCGTCTACTTTTCCC	-3'
59	consensus sequence	5'-	GYTTMRCTTCYRDGTTCG	-3'
60	consensus sequence	5'-	GCTTAACCTTCGTTGTTTCG	-3'
61	consensus sequence	5'-	GCTTAACCTTCTATGTTTCG	-3'
62	consensus sequence	5'-	GCTTAACCTTCTGTTGTTTCG	-3'
63	consensus sequence	5'-	GCTTAACCTTCCATGTTTCG	-3'
64	consensus sequence	5'-	GCTTAACCTCCGGGTTTCG	-3'
65	consensus sequence	5'-	GCTTAACCTCTAGGTTTCG	-3'
66	consensus sequence	5'-	GCTTAACCTCTGGGTTTCG	-3'
67	consensus sequence	5'-	GCTTAACCTCCAGGTTTCG	-3'
68	consensus sequence	5'-	GCTTAACCTCCGAGTTTCG	-3'
69	consensus sequence	5'-	GCTTAACCTCTAAGTTTCG	-3'
70	consensus sequence	5'-		

Table 2 (Cont.)

SEQ ID NO	Description		Sequence	
71	consensus sequence	5'-	GCTTAACTTCTGAGTTCG	-3'
72	consensus sequence	5'-	GCTTAACTTCCAAGTTCG	-3'
73	Lactobacillus	specific probe	TCGAGAAATAATTGAATAATATCTAG	-3'
74	Lactobacillus	specific probe	GAGGAAAGAAGTTCCTCTAT	-3'
75	Lactobacillus	specific probe	AACAGAGAAGATATTATCTAGTT	-3'
76	Lactobacillus	specific probe	TTGAGAGAACGAAGTTGCTCAGGCTTATGAAAAATAAGCAT	-3'
77	Lactobacillus	specific probe	TTTCGTTGGCCGGGTTTGGCCAAATGGATTTCAGGGTTCCTATGTGG	-3'
78	Lactobacillus	specific probe	GCCTTTTCGATGAATAACACTGGTTCCCGACAAACACAAAAACAACAATGA	-3'
79	Lactobacillus	specific probe	TTAGAAACCGGAGCATAAAGCGGCCCTGAG	-3'
80	Lactobacillus	specific probe	GCCTGATGCCGGGCTTTGGCCATTCGGTTCAGGTCCTTATGTGC	-3'
81	Lactobacillus	specific probe	CAAGTACGTTAAAGTTCAAGGCAGCAATTAACAAATAGTAGTAGTT	-3'
82	Lactobacillus	specific probe	AAAGAAATGAATATCCAGTTTTCAGAGCGCAACGTTCTCAGAAA	-3'
83	Lactobacillus	specific probe	AGGTGCAATGTTAGGCTTTTGAAATGAAATATTACTTATATATGCAGTT	-3'
84	Pediococcus	specific probe	GCCTGTAAGTGGATCGGAGAA	-3'
85	Pediococcus	specific probe	GCCGCGAAGTGGATCGGAGAA	-3'
86	Pediococcus	specific probe	GAGAGAATAAAATTTCTTTCACACGA	-3'
87	Pectinatus	specific probe	AAAATCATCGAAAAAAATGTTGGTCTGAGATTCTTCT	-3'
88	Pectinatus	specific probe	CACTCTGGTTGAAGGCAGGGAACG	-3'
89	Pectinatus	specific probe	GATTTTCATCAAAAAAGAGAAATGTTTGGTCAGAGATTTT	-3'
90	Pectinatus	specific probe	TATATACCGCTGAGGTGCTGAGGCACACTGAAG	-3'
91	Pectinatus	specific probe	AATTTTCATCTATAAAATGTTTGGTCTCCTGATTTCTTCT	-3'
92	Pectinatus	specific probe	AGATTAGTTCCTGGTTTACTTTATATATAGACACTAAGGTGCAGAAAAAG	-3'
93	Pectinatus	specific probe	AGGAAACGGCGCTTCGTAA	-3'
94	Selenomonas	specific probe	TAATAATCTAGAATGTTTCGATACAAATTTTCTTCTGTATAGTTTGTGAG	-3'
95	Zymophilus	specific probe	TGGACAT	-3'
96	Zymophilus	specific probe	GAGCGGAAAGCGGAAGGCAGCGAT	-3'
97	Megasphaera	specific probe	GAGCGGAAAGCTAAAGGCAGCGAT	-3'
		specific probe	AATCCTGAAACGAAATTCAGTGGTGTATGGCTGCAGGGA	-3'

Table 2 (Cont.)

SEQ ID NO	Description	Sequence	
98	Detection of all Lactobacillaceae relevant to brewing for differentiation from other bacteria relevant to brewing	TATGGAAGTAAGACCCCTGA	5'- -3'
99		AGATGATCAGGTAGTAGGCT	5'- -3'
100		AGATGATCAGGTCGATAGGTT	5'- -3'
101		AGATGATCAGGTAGATAGGTT	5'- -3'
102		TACTAATCGGTCGAGGACTTAACCA	5'- -3'
103		ATACTAATCAGTCGAGGACTTAACCA	5'- -3'
104		GAAGCGGACTGGTACTAATAAGCCGAGAGCTT	5'- -3'
105	Pectinatus	CAGCGGACCAATACTAATAAATCGAGGGCTTA	5'- -3'
106	Selenomonas	AGCGGACCGATACATAATAGGTCGAGGGCTTGACTTAAA	5'- -3'
107	Zymophilus	GGAGCGGACCGGTACTAATAGACCGAGGACTT	5'- -3'
	Megasphaera		

Table 3

[illegible]

	SEQ ID NO 33	SEQ ID NO 34	SEQ ID NO 35	SEQ ID NO 36	SEQ ID NO 37	SEQ ID NO 38	SEQ ID NO 39	SEQ ID NO 40-45
<i>Lactobacillus brevis</i>	-	-	-	-	-	-	-	+
<i>Lactobacillus lindneri</i>	-	-	-	-	-	-	-	+
<i>Lactobacillus casei</i>	-	-	-	-	-	-	-	+
<i>Lactobacillus paracasei paracasei</i>	-	-	-	-	-	-	-	+
<i>Lactobacillus coryniformis coryniformis</i>	-	-	-	-	-	-	-	+
<i>Lactobacillus coryniformis torquens</i>	-	-	-	-	-	-	-	+
<i>Lactobacillus curvatus</i>	-	-	-	-	-	-	-	+
<i>Pediococcus damnosus</i>	-	-	+	-	-	-	-	+
<i>Pediococcus inopinatus</i>	-	-	+	-	-	-	-	+
<i>Pectinatus cerevisiiphilus</i>	-	-	-	+	-	-	-	+
<i>Pectinatus frisingensis</i>	-	-	-	+	-	-	-	+
<i>Pectinatus sp. DSM 20462</i>	-	-	-	+	-	-	-	+
<i>Megasphaera cerevisiae</i>	-	-	-	-	+	-	-	+
<i>Selenomonas lacticifex</i>	-	-	-	-	-	+	-	+
<i>Zymophilus raffinosivorans</i>	+	-	-	-	-	-	+	+
<i>Zymophilus paucivorans</i>	-	+	-	-	-	-	+	+